



SCAN ME

Hybrid course

CRISPR Screening and Analysis

Wellcome Genome Campus, UK

Online learning dates: 28 July - 1 August | In-person training dates: 21-26 September

This two-phase course blends online learning with in-person, hands-on training to help you learn how to apply CRISPR screening and analysis techniques to real-world scenarios.

Experts from the Wellcome Sanger Institute and Open Targets will teach you everything you need to know about pooled and arrayed CRISPR screens.

The course will focus on three therapeutic areas where pooled and arrayed screens are commonly used:

- **Oncology**
- **Neurodegeneration**
- **Immunology**

Apply by 3 June to be considered for a place!



CRISPR - One of the biggest scientific stories
of the century

Gene Editing workflows and integrated concepts

Kalpana Surendranath MSc PhD PGCHE SFHEA FRSB

Associate professor in genome engineering

University of Westminster

Director - Gene Editors of the Future

www.geneeditorsofthefuture.uk



UNIVERSITY OF
WESTMINSTER




Learning objectives

At the end of the session, you are expected have an appreciation of:

- The origins
- Fundamentals of **CRISPR** /Cas9 tools
- Workflow of **CRISPR** in generating *cellular models of disease*-
Guide Design and Clonal isolation
- Benchling and vendor tools

CRISPR Timeline



<https://www.nature.com/articles/s41392-023-01309-7>

The molecular arms race

Altering the structure of their receptors

1 Adsorption block

2 Injection block

Bacteria can clog the phage entryway and block DNA injection

3 Replication block

Methylated Vs unmethylated DNA


4 RM systems

5 Cell suicide


abortive infection

6 CRISPR

Adaptive immunity





Bacterial adaptive immunity



CRISPR-conferred phage resistance relies on incorporation of short phage-derived sequences within specific loci of the bacterial genome. In effect, the bacterial cell becomes immune to phages which carry homologous sequences

CRISPR Locus




```

TATAAATCAGTAAGTTACGAGGCCTCGAAAAAGAGGGTTTCTGGCAGGAAAACTCGGTATTCCTTT
CCTTCAAATGGTTATAGGTTTCAGGGCTAGTTCAGTCCGTATAGGCAGCTAAGAAA GCAGCGATCAAC
GGCCATATCTACGAGCTGGA GTTCACTGCCGTATAGGCAGCTAAGAAATGACCCGCCGATGGTGCTGG
GCCTGCAGGA GTTCACTGCCGTATAGGCAGCTAAGAAATGCCAGCGGGCGGTATGCGCTGCGGAGCT
TC GTTCACTGCCGTATAGGCAGCTAAGAAATCGATCAGCTTCGCGGCGGGCGGTGATTCTTCACTG
CCACATAGGTCGTCAAGAAA CGGCCAGCCTGCGCCGTTTCATTGCCGTGTAGTCCCGTAGGGCGAATGCCG
  
```


- Prokaryotes generally have a combination of different Cas genes, but Cas1 and Cas2 are conserved.
- The Cas genes are located several hundred basepairs upstream of the CRISPR array
- The number of CRISPR repeat/spacer groups varies by species, but each repeat is always followed by a unique spacer

Programmable Nucleases




Meganuclease

20-40 bp/Enzyme



ZFN

3 bp/Finger




TALEN

1 bp/Module

High cost

difficult



CRISPR-Cas


1 bp/Base

As of July 28, 2025

MY CUSTOM FILTERS 

58,592 results

RESULTS BY YEAR



Did you mean ***crisp*** (4,329 results)



Genome modification by **CR**

1

Ma Y, Zhang L, Huang X.

Cite


FEBS J. 2014 Dec;281(23):5186-93. doi: 10.1002/febs.23400. PMID: 25315507

Free article.



Clustered regularly interspaced short palindromic repeats

<https://pubmed.ncbi.nlm.nih.gov/?term=CRISPR>


CRISPR-Cas complex basics



Guide RNAs



CRISPR / Cas9 complex



Nishimasu H, Ran FA, Hsu PD, et al. Crystal structure of Cas9 in complex with guide RNA and target DNA. *Cell*. 2014;156(5):935-949

Guide RNAs

By gRNA Type :




2023




2031

By Application :



By End-User:




Cas9 Variants

SpCas9 variant	Mutations (relative to SpCas9)	PAM sequence
D1135E variant	D1135E	NGG
VQR variant	D1135V, R1335Q and T1337R	NGAN or NGNG
EQR variant	D1135E, R1335Q and T1337R	NGAG
VRER variant	D1135V, G1218R, R1335E and T1337R	NGCG


Cas9 Orthologs

Cas9 species	PAM sequence ')
<i>Streptococcus pyogenes</i> (Sp)	3' NGG
<i>Staphylococcus aureus</i> (Sa)	NGRRT or NGRRN
<i>Neisseria meningitidis</i> (Nm or Nme)	NNNGATT
<i>Campylobacter jejuni</i> (Cj)	NNNNRYAC
<i>Streptococcus thermophilus</i> (St)	NNAGAAW
<i>Treponema denticola</i> (Td)	NAAAAC
~20 additional Cas9 species	PAM sequence may not be characterized


Key activities of Cas9



PAM occurrence on the human genome



B



Chr	NGG		NNAGAAW	
	median	mean	median	mean
1	7	12.8	67	115.8
2	8	12.7	64	100.8
3	8	13.0	63	98.5
4	9	14.0	61	94.5
5	8	13.1	63	97.9
6	8	13.1	63	98.5
7	8	12.4	64	102.9
8	8	12.8	64	100.9
9	7	13.9	65	120.5
10	7	12.1	66	107.0
11	7	12.0	65	105.8
12	8	12.4	65	103.5
13	8	13.6	62	94.6
14	8	12.0	65	101.5
15	7	11.5	68	107.7
16	7	11.7	74	136.8
17	6	10.3	76	127.9
18	8	13.4	63	101.8
19	6	9.4	82	145.4
20	7	11.1	72	121.8
21	7	13.4	64	111.4
22	6	9.2	85	140.3
X	8	13.2	63	99.0
Y	8	29.2	62	223.7



Design CRISPR guides: Guide parameters

Design type

- ☒ Single guide
Wild-type Cas9, single gRNA (higher efficiency)
- ☐ Paired guides
Double Cas9 nickase, two gRNAs (lower off-target effects)
- ☐ Guides for "base editing" (Komor et al., 2016)
C > T (or G > A) substitution, no dsDNA breaks

Guide length: 20

Genome: GRCh38 (hg38, Homo sapiens)

Don't see the genome you're looking for? We may be able to import it, just ask.

PAM: NGG (SpCas9, 3' side)

Custom PAM

- NGG (SpCas9, 3' side)
- NAG (SpCas9, 3' side)
- NG (SpCas9 NG, 3' side)
- NNNGATT (NmeCas9, 3' side)
- NNAGAAW (StCas9, 3' side)
- NAAAAC (TdCas9, 3' side)
- NNRR (SaCas9, 3' side)
- NNRRRT (SaCas9, 3' side)

☒ Save these as my default

☒ Explore Molecular

BASES 845

Finish

Cas9 proteins and PAMS

a


Cas9	PAM	Type
<i>Streptococcus pyogenes</i>	NGG	II-A
<i>Streptococcus mutans</i>	NGG	II-A
<i>Listeria innocua</i>	NGG	II-A
<i>Streptococcus thermophilus</i> A (CRISPR3)	NGGNG	II-A
<i>Francisella novicida</i>	NG	II-B
<i>Lactobacillus buchneri</i>	NAAAAN	II-A
<i>Treponema denticola</i>	NAAAAN	II-A
<i>Streptococcus thermophilus</i> B (CRISPR1)	NNAAAAW	II-A
<i>Campylobacter jejuni</i>	NNNNACA	II-C
<i>Pasteurella multocida</i>	GNNNCNNA	II-C
<i>Neisseria meningitidis</i>	NNNNGATT	II-C

b

PAM

<i>S.pyogenes</i>	1320	APAAFKYFD - - - - TT	IDRKR	-YTSTKEV - LDATL	IHQSI	ITG - -	LYETRIDL	SQLGGD - - - -	1368
<i>L.innocua</i>	1290	APASFKFFE - - - - TT	IERKR	-YNNLKEL - LNSTI	IYQSI	ITG - -	LYESRKRLDD - - - - -	1334	
<i>S.mutans</i>	1297	APATFKFFD - - - - KN	IDRKR	-YTSTTEI - LNATL	IHQSI	ITG - -	LYETRIDL	NKLGGD - - - -	1345
<i>S.thermophilus A</i>	1339	SAADFEEFLG - - - - VK	IPRYRDY	TPSSLL - KDATL	IHQSV	ITG - -	LYETRIDL	AKLGEG - - - -	1388
<i>F.novicida</i>	1569	DDSKINYFMN - - - - HS	LLKSRY	PDQVLEILKQSTI	IEFESSG	- -	FNKTIKEM	LGMKLAGI - -	1622
<i>L.buchneri</i>	1320	ANPTFGNLKD	IGITTPFGQL	QQPNGILLS - DEAKI	RYQSP	TG - -	LFERTVSL	KDL - - - - -	1371
<i>T.denticola</i>	1344	ATRNVS	DLQHIGGSKYSGVAK	IGNKISSL - DNCIL	IYQSI	ITG - -	IFEKRIDL	LLKV - - - - -	1395
<i>S.thermophilus B</i>	1078	NVANSGQCKKGLGKSN	ISIYKVRTDVL - -	GNQHI	IKNEGDK	- - - -	PKLDF - - - - -	1121	
<i>C.jejuni</i>	951	- - - - -	SIGIQNLKVFEKYIVS - - - -	ALGEVTKA	- - - -	EFRQREDFKK - - - - -	984		
<i>P.multocida</i>	1000	RATGNI	SLKEHDGEISK	GKDGVRVGVKL - - -	ALSFEKYQVD - -	ELGKNRQICRP	QQRQPVR	1056	
<i>N.meningitidis</i>	1025	RGTGNI	NIIRIHDLDHK	IGKNGILEGIGV - -	KTALSFQKYQID - -	ELGKEIRPCRL	KKRPPVR	1082	

Genome Diversity



Mistletoe (*Viscum album*): ~ 90 billion base pairs of DNA (bp), 10 chromosome pairs, 39,000 genes

Axolotl (*Ambystoma mexicanum*): ~ 32 billion base pairs, 14 chromosome pairs

Human (*Homo sapiens*): ~ 3.2 billion base pairs, 23 chromosome pairs, 21,000 genes




Horseshoe bat (*Rhinolophus ferrumequinum*): ~ 2.1 billion base pairs, 19,000 genes.

Chicken (*Gallus gallus*): ~ 1.2 billion base-pairs, 39 chromosomes (most of which are 'microchromosomes') 20,000 genes

The *Caenorhabditis elegans* worm: ~ 100 million base pairs, 6 chromosome pairs, 20,500 genes

Bacteria *Mycoplasma genitalium*: ~ 580,000 base pairs, 1 chromosome, 470 genes

The human genome 1997- 2003



APPROVED Clinical Investigation
Committee per: OK Consent Form, Page
Effective date: 3/24/97
Termination date: 1/21/98

Roswell Park Cancer Institute

Informed Consent from Blood Donor Chromosome (BAC) Libraries.

You are being asked to provide samples of blood and cell lines, to be used in research. This form so that you can decide whether or not you want

What we may do with your blood samples, w
To increase the protection of your identity as DNA approximately 10 different people for each donor means that you will not have any knowledge if the selection of the particular donor will be done your blood to prepare two kinds of blood products to use your blood for further research, then we components, cut the DNA in many fragments, and bacteria to make "Bacterial Artificial Chromosomes" or clones, of these BACs. The collection of all these as a "BAC library". The library we plan to make genetic material, distributed over many clones from the white blood cells from your blood and treat stay alive outside your body) and their ability to (cells multiplying outside your body, generating single white blood cell are designated as a "cell" prepared from your blood.

At present, BAC libraries, like the one we will use One is to find out where on the chromosomes information mapping. The second is to determine the sequence to finding out how DNA words that make up get other types of research in the future. [Chromosome cell; these structures contain an organized arrangement "blue-print". The basic segments are also designed information related to a particular hereditary trait structure of the chemical known as "DNA". The bars, the DNA consists of a string of beads ("base pairs") and the decoding of the hereditary message boils This process is known as "DNA Sequencing". This is the main reason why we prepare BAC libraries

protocol. [redacted] will attach a unique numerical label to your blood sample(s), but [redacted] will keep a list of names and corresponding numbers. The blood samples will then be transferred to Dr. De Jong and he will at random [blindly] select one of the blood samples from the different donors for further processing. He will not maintain the information on the original number of the sample, but will attach a new label to the tube (e.g. donor "n"). Dr. De Jong will not provide any feedback to [redacted] on his selection of a particular blood sample, to ensure that nobody will know that you have donated your blood for the preparation of a BAC library.

You will need to deal only with [redacted] when you have any questions about this research and this Informed-Consent form. The people who will take your blood sample will not know your identity. [redacted] will see that you are paid the agreed upon nominal amount only for your time and effort; it will be taken out of a separate cash fund so that you will not get checks directly from the Roswell Park Cancer Institute. Thus, your identity cannot be revealed through a financial auditing trail.

Since every person's DNA is unique, it will someday be possible that someone could find out who you are just from knowing your DNA sequence. At that time, the ease of identification would increase as more of a particular person's DNA sequence were known. To minimize the chance that someone would be able to find out who you are just from looking at the sequence in the public databases, several other DNA libraries are being created using DNA from other people. This means that even when the whole human DNA sequence is known, it will have been derived from DNA from a large number of people. If we use the blood you donate to generate a BAC library, then we expect that no more than 10% of the eventual DNA sequence will have been obtained from your DNA. However, more likely (approximately 90% probability), we will not have used your blood for the preparation of a BAC library. As a result of the precautions we have taken, the chance that someone could find out who you are just from knowing the DNA sequence, would be very low. In fact, no single person in our institute or elsewhere will have sufficient information to connect your identity with the BAC library. The principal investigator (Dr. De Jong) and his personnel will not have any access to the identity of the donors. [redacted]

Donor initials: _____


Date: _____

Witness initials: _____

Date: _____

RP11

Human Pangenome Reference Consortium



Freeze date	GENCODE release	Reference release?	Release date	Genome assembly version	UCSC version	Notes
03.2020	35	Y	08.2020	GRCh38	-	re-merge with new Havana annotation, updated Ensembl gene set
11.2019	34	N	04.2020	GRCh38	34	re-merge with new Havana annotation, updated Ensembl gene set
08.2019	33	N	01.2020	GRCh38	33	re-merge with new Havana annotation, updated Ensembl gene set
05.2019	32	N	09.2019	GRCh38	32 (current for GRCh38)	re-merge with new Havana annotation, updated Ensembl gene set
02.2019	31	N	06.2019	GRCh38	31	re-merge with new Havana annotation, updated Ensembl gene set
11.2018	30	N	04.2019	GRCh38	30	re-merge with new Havana annotation, updated Ensembl gene set
05.2018	29	N	10.2018	GRCh38	29	re-merge with new Havana annotation, updated Ensembl gene set

Human Pangenome Reference Consortium

Funded by: NIH NHGRI (National Human Genome Research Institute)




Collaboration Map




University of California, Santa Cruz (UCSC)
The Rockefeller University
The Broad Institute
Washington University in St. Louis
Baylor College of Medicine

Cancer Drug Discovery



Project **Achilles**



Seminal articles in CRISPR/Cas9



Published as: *Science*. 2012 August 17; 337(6096): 816–821.

A programmable dual RNA-guided DNA endonuclease in adaptive bacterial immunity

Martin Jinek^{1,2}, Krzysztof Chylinski^{3,4}, Ines Fonfara⁴, Michael Hauer^{2,5}, Jennifer A. Doudna^{1,2,6,7,*}, and Emmanuelle Charpentier^{4,*}

¹Howard Hughes Medical Institute, University of California, Berkeley, California 94720, USA. ²Department of Molecular and Cell Biology, University of California, Berkeley, California 94720, USA. ³Max F. Perutz Laboratories, University of Vienna, A-1030 Vienna, Austria. ⁴The Laboratory for Molecular Infection Medicine Sweden (MIMS), Umeå Centre for Microbial Research (UCMR), Department of Molecular Biology, Umeå University, S-90187 Umeå, Sweden. ⁵Present address: Friedrich Miescher Institute for Biomedical Research, 4058 Basel, Switzerland. ⁶Department of Chemistry, University of California, Berkeley, California 94720, USA. ⁷Physical Biosciences Division, Lawrence Berkeley National Laboratory, Berkeley, California 94720, USA.

* These authors contributed equally to this work.

Abstract

CRISPR/Cas systems provide bacteria and archaea with adaptive immunity against viruses and plasmids by using crRNAs to guide the silencing of invading nucleic acids. We show here that in a subset of these systems, the mature crRNA base-paired to *trans*-activating tracrRNA forms a two-RNA structure that directs the CRISPR-associated protein Cas9 to introduce double-stranded (ds) breaks in target DNA. At sites complementary to the crRNA-guide sequence, the Cas9 HNH nuclease domain cleaves the complementary strand while the Cas9 RuvC-like domain cleaves the non-complementary strand. The dual-tracrRNA:crRNA, when engineered as a single RNA chimera, also directs sequence-specific Cas9 dsDNA cleavage. Our study reveals a family of endonucleases that use dual RNA for site-specific DNA cleavage and highlights the structural

PROTOCOL

Genome engineering using the CRISPR-Cas9 system

F Ann Ran^{1–5,8}, Patrick D Hsu^{1–5,8}, Jason Wright¹, Vineeta Agarwal^{1,6,7}, David A Scott^{1–4} & Feng Zhang^{1–4}

¹Broad Institute of Massachusetts Institute of Technology (MIT) and Harvard, Cambridge, Massachusetts, USA. ²McGovern Institute for Brain Research, Cambridge, Massachusetts, USA. ³Department of Brain and Cognitive Sciences, MIT, Cambridge, Massachusetts, USA. ⁴Department of Biological Engineering, MIT, Cambridge, Massachusetts, USA. ⁵Department of Molecular and Cellular Biology, Harvard University, Cambridge, Massachusetts, USA. ⁶Program in Biophysics, Harvard University, MIT, Cambridge, Massachusetts, USA. ⁷Harvard-MIT Division of Health Sciences and Technology, MIT, Cambridge, Massachusetts, USA. ⁸These authors contributed equally to this work. Correspondence should be addressed to F.Z. (zhang@broadinstitute.org).

Published online 24 October 2013; doi:10.1038/nprot.2013.143

Targeted nucleases are powerful tools for mediating genome alteration with high precision. The RNA-guided Cas9 nuclease from the microbial clustered regularly interspaced short palindromic repeats (CRISPR) adaptive immune system can be used to facilitate efficient genome engineering in eukaryotic cells by simply specifying a 20-nt targeting sequence within its guide RNA. Here we describe a set of tools for Cas9-mediated genome editing via nonhomologous end joining (NHEJ) or homology-directed repair (HDR) in mammalian cells, as well as generation of modified cell lines for downstream functional studies. To minimize off-target cleavage, we further describe a double-nicking strategy using the Cas9 nickase mutant with paired guide RNAs. This protocol provides experimentally derived guidelines for the selection of target sites, evaluation of cleavage efficiency and analysis of off-target activity. Beginning with target design, gene modifications can be achieved within as little as 1–2 weeks, and modified clonal cell lines can be derived within 2–3 weeks.

INTRODUCTION

The ability to engineer biological systems and organisms holds enormous potential for applications across basic science, medicine and biotechnology. Programmable sequence-specific endonucleases that facilitate precise editing of endogenous genomic loci are now enabling systematic interrogation of genetic elements and causal genetic variations^{1–2} in a broad range of species, including those that have not previously been genetically tractable^{3–6}. A number of genome editing technologies have emerged in recent years, including zinc-finger nucleases (ZFNs)^{7–10}, transcription activator-like effector nucleases (TALENs)^{10–17} and the RNA-guided CRISPR-Cas nuclease system^{18–23}. The first two technologies use a strategy of tethering endonuclease catalytic domains to modular DNA-binding proteins for inducing targeted DNA double-stranded breaks (DSBs) at specific genomic loci. By contrast, Cas9 is a nuclease guided by small RNAs through Watson-Crick base pairing with target DNA^{24–28} (Fig. 1), representing a system that is markedly easier to design, highly specific, efficient and well-suited for high-throughput and multiplexed gene editing for a variety of cell types and organisms.

Precise genome editing using engineered nucleases

Similarly to ZFNs and TALENs, Cas9 promotes genome editing by stimulating a DSB at a target genomic locus^{29,30}. Upon cleavage by Cas9, the target locus typically undergoes one of two major pathways for DNA damage repair (Fig. 2): the error-prone NHEJ or the high-fidelity HDR pathway, both of which can be used to


frequencies than NHEJ), it can be leveraged to generate precise, defined modifications at a target locus in the presence of an exogenously introduced repair template. The repair template can either be in the form of conventional double-stranded DNA targeting constructs with homology arms flanking the insertion sequence, or single-stranded DNA oligonucleotides (ssODNs). The latter provides an effective and simple method for making small edits in the genome, such as the introduction of single-nucleotide mutations for probing causal genetic variations³¹. Unlike NHEJ, HDR is generally active only in dividing cells, and its efficiency can vary widely depending on the cell type and state, as well as the genomic locus and repair template³².

Cas9: an RNA-guided nuclease for genome editing

CRISPR-Cas is a microbial adaptive immune system that uses RNA-guided nucleases to cleave foreign genetic elements^{18–21,26}. Three types (I–III) of CRISPR systems have been identified across a wide range of bacterial and archaeal hosts, wherein each system comprises a cluster of CRISPR-associated (Cas) genes, noncoding RNAs and a distinctive array of repetitive elements (direct repeats). These repeats are interspaced by short variable sequences²⁰ derived from exogenous DNA targets known as protospacers, and together they constitute the CRISPR RNA (crRNA) array. Within the DNA target, each protospacer is always associated with a protospacer adjacent motif (PAM), which can vary depending on the specific CRISPR system^{34–36}.

The Type II CRISPR system is one of the best characterized^{25–28,37,38}

DNA damage




Nature Reviews | **Genetics**




Free radical species are generated either as a by-product of metabolism or through exposure to exogenous physical agents, such as ionizing radiation, which can induce the formation of double-strand breaks.


DNA repair pathways and Gene Editing




Double Strand Break damage




Homologous recombination Vs Non-Homologous end joining



Roadmap of CRISPR tools



Vendor options/tools



QUICK ORDER KALPANA SURENDRANATH



Search for Products, SKUs, Part #...

Search

0 MY CART

Products & Services

Explore Solutions

Support & Education

Tools

Company

Contact Us

Products › Oligos & Reagents › Custom DNA Oligos › DNA Oligonucleotides

Custom DNA Oligos

Consistently generate reliable data from high-quality oligos

All single-stranded and duplexed DNA sequences are produced with high coupling efficiencies, resulting in high-quality DNA products. Our specialized platforms allow us to deliver primers for PCR, dual-labelled probes for qPCR, indexed adapters and fusion primers for sequencing, and a variety of custom products.




Live Chat

Vendor options/tools

Alt-R gRNA OPTIONS FOR S.p. Cas9 & A.s. Cas12a


Guide RNAs		Cas9 guide RNAs			Cas12a guide RNAs	
Structure	Alt-R 2-part 	Alt-R 2-part XT 	Alt-R sgRNA 	Alt-R Cas12a crRNA 		
gRNA format	Alt-R CRISPR-Cas9 crRNA & tracrRNA	Alt-R CRISPR-Cas9 crRNA XT & tracrRNA	Alt-R CRISPR-Cas9 sgRNA	Alt-R CRISPR-Cas12a crRNA		
Components	crRNA tracrRNA	crRNA XT tracrRNA	sgRNA	crRNA		
Size (nt)	36 67	36 67	100	40–44 (41 nt recommended)		
Annealing required	Yes	Yes	No	No		
Stability	++	+++	++++	+++		
Applications	• Cas9-expressing cells • RNP in most cell types	• Co-delivery with Cas9 plasmid/Cas9 mRNA • RNP under difficult experimental conditions (e.g., high nuclease environments)		• KO/KI, RNP in most cell types • Cas12a-expressing cells		

Quantity	Product	Catalog #	Price
0	Alt-R™ S.p. Cas9 Nuclease V3, 500 µg	1081059	£653.00 GBP



Two-step enrichment of Cas9

Workflow for Knock out generation



Gene Editing Workflow- Stage 1

1. Design and selection of **guide sequences** (online algorithm)
2. Synthesis of single stranded guide oligos
3. Clone into **CRISPR/Cas9 expression vector**
4. Sanger Sequencing
5. Plasmids purification

Stage1


sgRNA synthesis and RNP complex formation

6. Transfect cells
7. Selection e.g. **Puromycin**
8. Clonal Isolation and expansion

Stage2 ~ 2 months

Stage 3

9. Functional characterization and phenotype and genotype analysis



Plasmids tools

Search

We found 23,334 results for: crispr ✕

Select a Category to Narrow Your Search

Catalog

- 22,561 in Plasmids
- 16 in Viral Preps
- 270 in Pooled Libraries
- 27 in Plasmid Kits
- 12 in Bacterial
- [Help Center](#)

Early days...

<http://crispr.mit.edu>

Haeussler, M., Schönig, K., Eckert, H. *et al.* Evaluation of off-target and on-target scoring algorithms and integration into the guide RNA selection tool CRISPOR. *Genome Biol* **17**, 148 (2016).

Guide Design Webtools

Algorithm	Editing	URL
Cas-Designer	CRISPR-Cas mediated editing	http://www.rgenome.net/cas-designer/
Cas-OFFinder		http://www.rgenome.net/cas-offfinder/
CHOPCHOP		https://chopchop.cbu.uib.no/
CRISPR-P 2.0		http://crispr.hzau.edu.cn/CRISPR2/
CRISPOR		http://crispor.tefor.net/
CRISPRdirect		https://crispr.dbcls.jp/
E-CRISP		http://www.e-crisp.org/
CCTop		http://crispr.cos.uni-heidelberg.de/
DeepCpf1		http://deepcrispr.info/
DeepSpCas9		http://deepcrispr.info/DeepSpCas9
DeepHF		http://www.DeepHF.com/
CINDEL		http://big.hanyang.ac.kr/cindel
SSC		http://cistrome.org/SSC/
Prediction web tool		https://crispr.ml
CRISPRscan		https://www.crisprscan.org
Microhomology-Predictor		http://www.rgenome.net/michcalculator
FORECasT		https://partslab.sanger.ac.uk/FORECasT
inDelphi		https://indelphi.giffordlab.mit.edu
TIDE		http://shinyapps.datacurators.nl
Cas-Analyzer		http://www.rgenome.net/cas-analyzer/
CRISPResso2		https://crispresso.pinellolab.partners.org/
CRISPRAnalyzerR		http://crispr-analyzer.dkfz.de/
CRISPR-Sub		http://www.rgenome.net/crispr-sub
BE-Designer	Base editing	http://www.rgenome.net/be-designer/
DeepBaseEditor		http://deepcrispr.info/DeepBaseEditor
BE-Hive		https://www.crisprbehive.design/
BE-Analyzer	Prime editing	http://www.rgenome.net/be-analyzer/
PE-Designer		http://www.rgenome.net/pe-designer
PrimeDesign		https://drugthatgene.pinellolab.partners.org/
PE-Analyzer		http://www.rgenome.net/pe-analyzer/
DeepPE		http://deepcrispr.info/DeepPE/
pegFinder		http://pegfinder.sidichenlab.org

Guide Verification tools - IDT

☐ Checker analysis: **Sequence1**

Custom checker Alt-R CRISPR-Cas9 gRNA

Sequence	On-target score	Off-target score
TGTCTCGCGAGCTCAGAGCG	22	89

[Hide off-target details -](#) | [Show related products](#)

+ ADD TO DESIGN SET


CREATE NEW DESIGN

Off-target results Demonstrates off-target hits to corresponding gRNA Design. [Export to Excel](#) | [Export to CSV](#)

Sequence	PAM	Score	#MM	Gene	Locus
TGTCTCGCGAGCTCAGAGCG	GGG	N/A		ZFP36L1	chr14:+68790299
TGTCACAGCAGCTCAGAGCG	CAG	21	4		chr18:-251793
TGTTCATGAGCTCAGAGTG	TGG	31	4		chr1:-60293221
TGTCTCGC-AGCTCAGAGCA	TGG	36	2		chr17:-78526399
AGTCTCGCGAGCCGGAGCG	TGG	37	3	NPRL3	chr16:-138651
TGTCTCG-GAGCTCAGGGCG	CAG	38	2	ARF3	chr12:+48937327
CTGTCGCTGGCTCAGAGCG	GGG	43	4		chr2:-240820638
TTTCTCAAAGAGCTCAGAGCT	AAG	46	4		chr5:-39232040
TGTCTCTCGAACTCGGAGCT	GGG	46	4	ZC3HC1	chr7:-130023680

https://eu.idtdna.com/site/order/designtool/index/CRISPR_SEQUENCE

Guide Verification tools - Synthego




OFF TARGET SEQUENCES						Filters
OFF TARGET SITE ?	MISMATCHES ?	▲ CHROMOSOME ?	CUT SITE ?	PAM ?	GENE ?	
AGTCTCGCGAGCCCGGAGCG	3	16	138,655	TGG	NPRL3	
TGTCTCGTGAGGTTAGAGTG	4	11	34,667,128	AGG		
TGTCTGTCAAGCTCAGAGCA	4	7	44,384,737	TGG	NUDCD3	
TGTCTCGGAAACTCAGAGCC	4	5	112,861,308	GGG	CTC-554D6.1	
TGTCTCCCCAGCTCACAGCC	4	CHR_HG30_PATCH	179,903,182	TGG		
TGTCTCCAAGCTCAGGGCC	4	10	93,070,733	TGG		
TGTCTCTGGAGCTCAGCGCA	4	4	111,753,901	TGG		
AGGCACGCGAGTTCAGAGCG	4	GL000205.2	56,922	TGG		

Construction of the guide expression plasmids


Stage 1



Bacterial Colony PCR



Plasmids sequence verification





BbsI overhang

gRNA 3 insert

TTTATATATCTTGTGNAAGGACGAAACACCGGTCTCGCGAGCTCAGAGCGGGTTTTTA
GAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCA
CCGAGTCGGTGCTTTTTTTGTTTTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCC
GTTTTTAGCGCGTGCGCCAATTCTGCAGACAAATGGCTCTAGAGGTACCCGTTACATA
ACTTACGGTAAATGGCCCGCCTGGCTGACCGCCCAACGACCCCGCCCAATTGACGT
CAATAGTAACGCCAATAGGGACTTTCCATTGACGTCAATGGGTGGAGTATTTACGGTAA
ACTGCCCACTTGGCAGTACATCAAGTGTATCATATGCCAAGTACGCCCCCTATTGACG
TCAATGACGGTAAATGGCCCGCCTGGCATTGTGCCAGTACATGACCTTATGGGACTT
TCCTACTTGGCAGTACATCTACGTATTAGTCATCGCTATTACCATGGTCGAGGTGAGCC
CCACGTTCTGCTTCACTCTCCCCATCTCCCCCCCCCTCCCCACCCCAATTTTGTATTT
ATTTATTTTTTAATTATTTTGTGCAGCGATGGGGGCGGGGGGGGGGGGGGGGGGGTTTG


U6- Forward primer

Gene Editing workflow - Stage 2





Transfection
Selection
Clonal Isolation
Clonal expansion
Functional Characterization

Western blot



Sanger sequencing



Gene Editing Workflow- Stage 2& 3



Stage2
~ 2 months

- 6. Transfect cells
- 7. Selection e.g. **Puromycin**
- 8. Clonal Isolation and expansion





Stage 3

- 9. Functional characterization and phenotype and genotype analysis

Colony picking




Colony picking



Clonal selection – observe early!


HEK293



A clone of ~10cells





HEK293




A clone of ~100cells

Antibiotic selection – not too early!



PCR genotyping




➔ **Design genotyping primers during the guide RNA design**

<https://www.genewiz.com/en-gb/public/services/next-generation-sequencing/amplicon-sequencing-services>

<https://www.genewiz.com/en-gb/public/services/sanger-sequencing/sanger-ez>


PCR genotyping



← Amplicon expected
size 440 bp

ATCTCTGCTCCTATTGGCTGGATATTTTCGTATTCCCCGAGCTCCTAA
AAACG *gPCRFor*
*AACCAATAGGAAGAGCGGAC*AGCGATCTCTAACGCGCAAGCGCATA
TCCTTCTAGGTAGCGGGCAGTAGCCGCTTCAGGGAGGGACGAA
GAGACCCAGCAACCCACAGAGTTGAGAAATTTGACTGGCATTCAA
GCTGTCCAATCAATAGCTGCCGCTGAAGGGTGGGGCTGGATGGCG
TAAGCTACAGCTGAAGGAAGAACGTGAGCACGAGGCACTGAGGT
GATTGGCTGAAGGCACTTCCGTTGAGCATCTAGACGTTTCCTTGGC
TCTTCTGGCGCCAAAATGTCGTTCGTGGCAGGGGTATTTCGGCGG
CTGGACGA *GACAGTGGTGAACCGCATCGCGGCGGGGGAAGTTATC*
*CAG*CGGCCAGCTAATGCTATCAAAGAGATGATTGAGAACTGGTAC
GGAGGGAGTCGAGCCG *GGCTCACTTAAGGGCTACGA gPCRRev*
CTTAACGGGCCGCGTCACTCAATGGCGCGGACACGCCTCTTTGCC
CGGGCAGAGGCATGTACAGCGCATGCCACAACGGCGGAGGCCG
CCGGGTTCCCTGACGTGCCAGTCAGGCCTTCTCCTTTTCCGCAGA
CCGTGT

Gene Editing Verification



Verification of Gene Editing in Mixed pools and Monoclonal lines

in CRISPRESSO session by Khalid Akram

Gene Editing Verification



A rectangular sign with a thin, light-colored wooden frame is centered on a light gray background. The sign has a white background and features the text "TIME TO TAKE A BREAK" in a bold, black, italicized sans-serif font. The text is arranged in three lines: "TIME TO" on the top line, "TAKE A" on the middle line, and "BREAK" on the bottom line.

***TIME TO
TAKE A
BREAK***

Wellcome Sanger Institute Genome Editing (WGE) is a website that provides tools to aid with genome editing of human and mouse genomes

CRISPR Finder [Ensembl for Human](#) [Ensembl for Mouse](#)

The CRISPR Finder will show CRISPR sites (paired or single) in and around genes. You can ask the finder to score the pairs for potential off-target sites, and browse individual and paired CRISPR sites using the Genovise genome browser tool. We also provide the ability to find CRISPRs in genomic sequence or by gRNA.

Find CRISPRs in our genome browser:



Find CRISPRs by gene using our table:

Ensembl ID	Spencer	Status	Summary
ENSG00000271017	1	Complete	chr1:1,000,000-1,000,000

Find CRISPRs by 20bp gRNA:

Sequence: ATATGATGACATTAAGTCTCT

Species: ☐ Human (HOMO) ☒ Mouse (MUS) ☐ Rat (RAT)

Find CRISPRs

Find CRISPRs in genomic sequence:

Search Again

CRISPR ID: ENSG00000271017

chr1:1,000,000-1,000,000

No

Find off-targets by sequence:

Find Off-Targets

Sequence: ATATGATGACATTAAGTCTCT

chr1:1,000,000-1,000,000

No

If you use this site in your research, please cite:
WGE: A CRISPR database for genome engineering. Alex Hodgkins; Anna Farne; Sajith Perera; Tiago Grego; David J. Parry-Smith; William C. Skarnes; Vivek Iyer (Bioinformatics 2015)
doi:10.1093/bioinformatics/btv308

Copyright (c) 2019 Genome Research Limited Your use of this site indicates your agreement to the GNU AGPLv3 licence



CRISPick

Ranks and picks candidate CRISPRko/a/i sgRNA sequences to maximize on-target activity for target(s) provided.

GENETIC PERTURBATION PLATFORM

Last update: 2025-07-17 | [Changelog](#)

NEW (2024-09-24)! For the GRCh38 genome, CRISPick has added support for [gnomAD](#) variation filtering. This feature is enabled by default, and deprioritizes selection of sgRNAs that overlap sites of variation, e.g. SNPs and indels. See the "Advanced Options" section below for details.

Reference Genome

☐ Human GRCh38 (NCBI RefSeq vGCF_000001.405.40-RS_2024.08)

☐ Human GRCh38 (Ensembl v.114)

☐ Human GRCh37 (NCBI RefSeq vGCF_000001.405.25-RS_2024.09)

☐ Mouse GRCh38 (NCBI RefSeq v.108.20200622)

☐ Mouse GRCh38 (Ensembl v.102)

☐ Rat mRatBN7.2 (NCBI RefSeq vGCF_015227675.2-RS_2023.06)

☐ Rat mRatBN7.2 (Ensembl v.113)

Mechanism

☐ CRISPRko

☐ CRISPRa

☐ CRISPRI

Full-genome CRISPick downloads

How CRISPick works

FAQ



CRISPOR ([citation](#)) is a program that helps design, evaluate and clone guide sequences for the CRISPR/Cas9 system. [CRISPOR Manual](#)
See [full list of changes](#)

Step 1

Planning a lentiviral gene knockout screen? Use [CRISPOR Batch](#)

Sequence name (optional):

Enter a single genomic sequence, < 2300 bp, typically an exon

Clear Box - Reset to default

cttcctttgtccccaatctggggcgcgccggcgcccccctggcgccctaaggactcggcgccgggaagtggccaggcgggggcgacctcggtctacagcgcccggtattctcgcagctcaccatgGATGATGATATCGCCGCGCTCGTCGTGCAACGGCTCCGGCATGTGCAAGGCCGGCTTCGCGGGCGACGATGCCCGCCGGCCGTCTTCCCTCCATCGTGGGGCGCC

Step 2

Select a genome

Homo sapiens - Human - UCSC Feb. 2009 (GRCh37/hg19) + SNPs: 1000Genomes, ExaC

Note: pre-calculated exonic guides for this species are on the [UCSC Genome Browser](#). We have 1297 genomes, but not yours? Search [NCBI assembly](#) and send a GCF/_GCA_ ID to [CRISPOR support](#).

Step 3

Select a Protospacer Adjacent Motif (PAM)

20bp-NGG - Sp Cas9, SpCas9-HF1, eSpCas9 1.1

See [notes on enzymes](#) in the manual.

SUBMIT

Text case is preserved, e.g. you can mark ATGs with lowercase. Instead of a sequence, you can paste a chromosome range, e.g. chr1:11,130,540-11,130,751

Version 5.2 - [Documentation](#) - [Contact us](#) - [Downloads/local installation](#) - [Citation](#) - [License](#)



CRISPR Guide RNA Design Tool

The CRISPR Guide RNA design tool allows you to visualize, optimize, and annotate multiple gRNA sequences at a time. Get on and off-target scores in seconds to compare and optimize for higher activity and lower off-target effects. Save, organize and attach guides to relevant sequences and assemble gRNAs into plasmids — all within Benchling.

Sign up for free

Request a demo



Useful links

- <https://www.sanger.ac.uk/tool/wge-crispr-design-tool/>
- <https://www.opentargets.org/>
- <https://cellmodelpassports.sanger.ac.uk/>
- <https://depmap.org/portal/achilles/>
- <https://sites.broadinstitute.org/ccle/>
- <https://depmap.org/portal/>
- <https://depmap.sanger.ac.uk/>

All the best!

